

## Product Datasheet

### CPI-0610 carboxylic acid (orb1179526)

<b>Description</b>	CPI-0610 carboxylic acid
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{20}H_{15}ClN_2O_3$
<b>Purity</b>	> 98%
<b>MW</b>	366.8
<b>CAS Number</b>	[1380089-81-5]
<b>SMILES</b>	<chem>O=C(O)C[C@@H]1N=C(C2=CC=C(Cl)C=C2)C3=CC=CC=C3C4=C1O=C4C</chem>
<b>Expiration Date</b>	12 months from date of receipt.