

Product Datasheet

LCL521 dihydrochloride (orb1179551)

Description	LCL521 dihydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{31}H_{54}Cl_2N_4O_7$
Purity	> 98%
MW	665.69
CAS Number	[1226759-47-2]
SMILES	<chem>O=C(CN(C)C)O[C@H](C1=CC=C([N+][O-])=O)C=C1[C@@H](COC(CN(C)C)=O)NC(CCCCCCCCCCCC)=O.[H]Cl.[H]Cl</chem>
Expiration Date	12 months from date of receipt.