

Product Datasheet

cis-Pralsetinib (orb1179592)

Description	cis-Pralsetinib
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₇ H ₃₂ FN ₉ O ₂
Purity	> 98%
MW	533.6
CAS Number	[2097132-93-7]
SMILES	<chem>O=C([C@@]1(OC)CC[C@@H](C2=NC(NC3=NNC(C)=C3)=CC(C)=N2)CC1)N[C@H](C4=CC=C(N5N=CC(F)=C5)N=C4)C</chem>
Expiration Date	12 months from date of receipt.