

Product Datasheet

ENMD-2076 Tartrate (orb1179595)

Description	ENMD-2076 Tartrate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{31}N_7O_6$
Purity	> 98%
MW	525.56
CAS Number	[1291074-87-7]
SMILES	<chem>O[C@H](C(O)=O)[C@@H](C(O)=O)O.CC1=CC(NC2=NC(/C=C/C3=CC=C/C=C3)=NC(N4CCN(CC4)C)=C2)=NN1</chem>
Expiration Date	12 months from date of receipt.