

## Product Datasheet

### GDC-0575 dihydrochloride (orb1179606)

<b>Description</b>	GDC-0575 dihydrochloride
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{23}BrCl_2N_4O$
<b>Purity</b>	> 98%
<b>MW</b>	450.2
<b>CAS Number</b>	[1657014-42-0]
<b>SMILES</b>	<chem>O=C(C1CC1)NC2=CNC3=C2C(N4C[C@H](N)CCC4)=C(Br)C=C3.[H]Cl.[H]Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.