

Product Datasheet

(S)-(-)-Citronellal (orb1179627)

Description	(S)-(-)-Citronellal
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₁₀ H ₁₈ O
Purity	> 98%
MW	154.25
CAS Number	[5949-05-3]
SMILES	C/C(C)=CCC[C@H](C)CC=O
Expiration Date	12 months from date of receipt.