

## Product Datasheet

### (Rac)-PT2399 (orb1179706)

<b>Description</b>	(Rac)-PT2399
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{10}F_5NO_4S$
<b>Purity</b>	> 98%
<b>MW</b>	419.32
<b>CAS Number</b>	[1672662-07-5]
<b>SMILES</b>	<chem>N#CC1=CC(F)=CC(OC2=CC=C(S(=O)(C(F)F)=O)C3=C2CC(F)(F)C3O)=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.