

Product Datasheet

PhiKan 083 hydrochloride (orb1179709)

Description	PhiKan 083 hydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{16}H_{19}ClN_2$
Purity	> 98%
MW	274.79
CAS Number	[1050480-30-2]
SMILES	<chem>CNCC1=CC2=C(C=C1)N(CC)C3=C2C=CC=C3.[H]Cl</chem>
Expiration Date	12 months from date of receipt.