

Product Datasheet

Bay-K-8644 ((R)-(+)-) (orb1179830)

Description	Bay-K-8644 ((R)-(+)-)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{16}H_{15}F_3N_2O_4$
Purity	> 98%
MW	356.3
CAS Number	[98791-67-4]
SMILES	<chem>CC1=C([N+])([O-])=O)[C@H](C2=C(C(F)(F)F)C=CC=C2)C(C(OC)=O)=C(C)N1</chem>
Expiration Date	12 months from date of receipt.