

Product Datasheet

Centrinone-B (orb1179977)

Description	Centrinone-B
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₇ H ₂₇ F ₂ N ₇ O ₅ S ₂
Purity	> 98%
MW	631.67
CAS Number	[1798871-31-4]
SMILES	FC(C([N+])([O-])=O)=CC=C1)=C1CS(C2=CC(F)=C(SC3=NC(NC4=NNC(C)=C4)=C(OC)C(N5CCCCC5)=N3)C=C2)(=O)=O
Expiration Date	12 months from date of receipt.