

Product Datasheet

FPR Agonist 43 (orb1180050)

Description	FPR Agonist 43
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{20}H_{21}ClN_4O_2$
Purity	> 98%
MW	384.86
CAS Number	[903895-98-7]
SMILES	<chem>O=C(NC1=C(C(C)C)N(C)N(C2=CC=CC=C2)C1=O)NC3=CC=C(Cl)C=C3</chem>
Expiration Date	12 months from date of receipt.