

Product Datasheet

LXR β agonist-2 (orb1180218)

Description	LXR β agonist-2
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₃₂ H ₃₁ F ₆ N ₃ O ₇
Purity	> 98%
MW	683.59
CAS Number	[1949801-52-8]
SMILES	<chem>OC(C(F)(F)F)(C(F)(F)F)C1=CC(CCC)=C(OC2=CC(O)=C(C(CN3C(NC(C)(C4=NC=C(OC(C)C)C=C4)C3=O)=O)=O)C=C2)C=C1</chem>
Expiration Date	12 months from date of receipt.