

## Product Datasheet

### 5-Acetylsalicylic acid (orb1180756)

<b>Description</b>	5-Acetylsalicylic acid
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>9</sub> H <sub>8</sub> O <sub>4</sub>
<b>Purity</b>	> 98%
<b>MW</b>	180.16
<b>CAS Number</b>	[13110-96-8]
<b>SMILES</b>	<chem>O=C(O)C1=CC(C(C)=O)=CC=C1O</chem>
<b>Expiration Date</b>	12 months from date of receipt.