

Product Datasheet

Cyclophilin inhibitor 1 (orb1180788)

Description	Cyclophilin inhibitor 1
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₃₁ H ₃₉ N ₅ O ₇
Purity	> 98%
MW	593.67
CAS Number	[1509904-96-4]
SMILES	<chem>C[C@@H](NC([C@@H]1CCCN(C([C@H](C)NC([C@H](C(C)C)OC2=O)=O)=O)N1)=O)C3=CC=C4C(C=C/C=C/C52COCOC5)C=C4)=N3</chem>
Expiration Date	12 months from date of receipt.