

Product Datasheet

cIAP1 Ligand-Linker Conjugates 3 (orb1181199)

Description	cIAP1 Ligand-Linker Conjugates 3
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{39}H_{56}N_4O_{11}S$
Purity	> 98%
MW	788.95
CAS Number	[2222354-20-1]
SMILES	<chem>C[C@H](N(C(OC(C)(C)C)=O)C)C(N[C@@H](C1CCCCC1)C(N(CCC2)[C@@H]2)C3=NC(C4=CC=CC(OCCOCCOCCOCC(O)=O)=C4)=O)=CS3)=O)=O</chem>
Expiration Date	12 months from date of receipt.