

## Product Datasheet

### Wnt/ $\beta$ -catenin agonist 1 (orb1181411)

<b>Description</b>	Wnt/ $\beta$ -catenin agonist 1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>22</sub> H <sub>25</sub> N <sub>3</sub> O <sub>2</sub>
<b>Purity</b>	> 98%
<b>MW</b>	363.45
<b>CAS Number</b>	[2305372-67-0]
<b>SMILES</b>	<chem>CCN1CC2=CC(OCC3=CC(C)=NN3C)=C(O)C=C2C4=C1C(C)=CC=C4</chem>
<b>Expiration Date</b>	12 months from date of receipt.