

## Product Datasheet

# DL-alpha-Tocopherol methoxypolyethylene glycol succinate (orb1181440)

<b>Description</b>	DL-alpha-Tocopherol methoxypolyethylene glycol succinate
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$(C_2H_4O)_n C_{34}H_{56}O_5$
<b>Purity</b>	> 98%
<b>CAS Number</b>	[1309573-60-1]
<b>SMILES</b>	<chem>CC1=C(CCC(CCCC(C)CCCC(C)CCCC(C)C)(C)O2)C2=C(C)C(C)=C1OC(CCC(OCCOC)=O)=O.[n]</chem>
<b>Expiration Date</b>	12 months from date of receipt.