

## Product Datasheet

### Phortress (orb1181582)

<b>Description</b>	Phortress
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{20}H_{25}Cl_2FN_4OS$
<b>Purity</b>	> 98%
<b>MW</b>	459.4
<b>CAS Number</b>	[741241-36-1]
<b>SMILES</b>	<chem>NCCCC[C@H](N)C(NC1=CC=C(C2=NC(C=C(F)C=C3)=C3S2)C=C1C)=O.Cl.Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.