

## Product Datasheet

### Piperidolate hydrochloride (orb1181728)

<b>Description</b>	Piperidolate hydrochloride
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>21</sub> H <sub>26</sub> ClNO <sub>2</sub>
<b>Purity</b>	> 98%
<b>MW</b>	359.89
<b>CAS Number</b>	[129-77-1]
<b>SMILES</b>	<chem>O=C(OC1CN(CC)CCC1)C(C2=CC=CC=C2)C3=CC=CC=C3.[H]Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.