

Product Datasheet

Bivalirudin TFA (orb1181755)

Description	Bivalirudin TFA
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{98}H_{137}DN_{24}O_{33} \cdot C_2HF_3O_2$
Purity	> 98%
MW	2294.31
CAS Number	[1191386-55-6]
SMILES	<chem>O=C(O)C(F)(F)F.[{d-Phe}-PRPGGGGNGDFEEIPEEYL]</chem>
Expiration Date	12 months from date of receipt.