

Product Datasheet

Ibrutinib Racemate (orb1181772)

Description	Ibrutinib Racemate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{24}N_6O_2$
Purity	> 98%
MW	440.5
CAS Number	[936563-87-0]
SMILES	<chem>C=CC(N1CC(N2N=C(C3=CC=C(OC4=CC=CC=C4)C=C3)C5=C(N)N=CN=C52)CCC1)=O</chem>
Expiration Date	12 months from date of receipt.