

Product Datasheet

Haloperidol Decanoate (orb1181798)

Description	Haloperidol Decanoate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{31}H_{41}ClFNO_3$
Purity	> 98%
MW	530.12
CAS Number	[74050-97-8]
SMILES	<chem>CCCCCCCCC(OC1(C2=CC=C(Cl)C=C2)CCN(CCCC(C3=CC=C(F)C=C3)=O)CC1)=O</chem>
Expiration Date	12 months from date of receipt.