

Product Datasheet

Tartaric acid (orb1181809)

Description	Tartaric acid
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₄ H ₆ O ₆
Purity	> 98%
MW	150.08
CAS Number	[87-69-4]
SMILES	<chem>O=C(O)[C@H](O)[C@@H](O)C(O)=O</chem>
Expiration Date	12 months from date of receipt.