

Product Datasheet

Emedastine Difumarate (orb1181848)

Description	Emedastine Difumarate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{34}N_4O_9$
Purity	> 98%
MW	534.56
CAS Number	[87233-62-3]
SMILES	<chem>CCOCCN1C(=NC2=CC=CC=C12)N3CCN(CCC3)C.C(=C/C(=O)O)C(=O)O.C(=C/C(=O)O)C(=O)O</chem>
Expiration Date	12 months from date of receipt.