

Product Datasheet

L-Azetidine-2-carboxylic acid (orb1181860)

Description	L-Azetidine-2-carboxylic acid
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₄ H ₇ NO ₂
Purity	> 98%
MW	101.11
CAS Number	[2133-34-8]
SMILES	C1CN[C@@H]1C(=O)O
Expiration Date	12 months from date of receipt.