

Product Datasheet

Amodiaquine hydrochloride (orb1181862)

Description	Amodiaquine hydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{20}H_{24}Cl_3N_3O$
Purity	> 98%
MW	428.78
CAS Number	[69-44-3]
SMILES	<chem>CCN(CC1=CC(NC2=CC=NC3=CC(Cl)=CC=C23)=CC=C1O)CC.[H]Cl.[H]Cl</chem>
Expiration Date	12 months from date of receipt.