

Product Datasheet

Cefodizime sodium (orb1181923)

Description	Cefodizime sodium
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{20}H_{18}N_6Na_2O_7S_4$
Purity	> 98%
MW	628.63
CAS Number	[86329-79-5]
SMILES	<chem>O=C(C(N12)=C(CSC3=NC(C)=C(CC(O[Na])=O)S3)CS[C@]2([H])[C@H](NC(/C(C4=CSC(N)=N4)=NOC)=O)C1=O)O[Na]</chem>
Expiration Date	12 months from date of receipt.