

## Product Datasheet

### Piperonyl butoxide (orb1181931)

<b>Description</b>	Piperonyl butoxide
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>19</sub> H <sub>30</sub> O <sub>5</sub>
<b>Purity</b>	> 98%
<b>MW</b>	338.44
<b>CAS Number</b>	[51-03-6]
<b>SMILES</b>	<chem>CCCC1=C(COCCOCCOCCCC)C=C2OCOC2=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.