

Product Datasheet

TCS-OX2-29 HCl (orb1182029)

Description	TCS-OX2-29 HCl
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{32}ClN_3O_3$
Purity	> 98%
MW	433.98
CAS Number	[1610882-30-8]
SMILES	<chem>CC(C)(C)[C@H](NCC1=CC=NC=C1)C(N2CC3=C(C=C(OC)C(OC)=C3)CC2)=O.[H]Cl</chem>
Expiration Date	12 months from date of receipt.