

## Product Datasheet

### FAA1 agonist-1 (orb1182084)

<b>Description</b>	FAA1 agonist-1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>21</sub> H <sub>17</sub> ClO <sub>5</sub> S
<b>Purity</b>	> 98%
<b>MW</b>	416.87
<b>CAS Number</b>	[2102196-57-4]
<b>SMILES</b>	<chem>C1C=CC=CC=C1C2=CC(COC3=CC=C(S(CC(O)=O)(=O)=O)C=C3)=CC=C2</chem>
<b>Expiration Date</b>	12 months from date of receipt.