

Product Datasheet

AZD1208 HCl (orb1182121)

Description	AZD1208 HCl
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{21}H_{22}ClN_3O_2S$
Purity	> 98%
MW	415.94
CAS Number	[1621866-96-3]
SMILES	<chem>O=C(NC(=O)SC1=C/C2=C(N3C[C@H](N)CCC3)C(C4=CC=CC=C4)=CC=C2.[H]Cl</chem>
Expiration Date	12 months from date of receipt.