

Product Datasheet

Ribocil B (orb1182130)

Description	Ribocil B
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{22}N_6O_5$
Purity	> 98%
MW	382.48
CAS Number	[1825355-55-2]
SMILES	<chem>O=C1N=C([C@H]2CCCN(CC3=CN=C(NC)N=C3)C2)NC(C4=CC=CS4)=C1</chem>
Expiration Date	12 months from date of receipt.