

## Product Datasheet

### Fenoprofen Calcium hydrate (orb1182164)

<b>Description</b>	Fenoprofen Calcium hydrate
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{15}H_{15}Ca_0.5O_4$
<b>Purity</b>	> 98%
<b>MW</b>	279.32
<b>CAS Number</b>	[71720-56-4]
<b>SMILES</b>	<chem>CC(C([O-])=O)C1=CC(OC2=CC=CC=C2)=CC=C1.O.[0.5Ca2+]</chem>
<b>Expiration Date</b>	12 months from date of receipt.