

Product Datasheet

EIDD-1931 (orb1182165)

Description	EIDD-1931
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_9H_{13}N_3O_6$
Purity	> 98%
MW	259.22
CAS Number	[3258-02-4]
SMILES	<chem>O[C@@H]1[C@@H](CO)O[C@](N2C(N=C(NO)C=C2)=O)([H])[C@@H]1O</chem>
Expiration Date	12 months from date of receipt.