

Product Datasheet

SR9011 hydrochloride (orb1182218)

Description	SR9011 hydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{32}Cl_2N_4O_3S$
Purity	> 98%
MW	515.5
SMILES	<chem>CCCCCNC(=O)N1CCC(C1)CN(CC2=CC=C(C=C2)Cl)CC3=CC=C(S3)[N+](=O)[O-].Cl</chem>
Expiration Date	12 months from date of receipt.