

Product Datasheet

Batimastat sodium salt (orb1182406)

Description	Batimastat sodium salt
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{30}N_3NaO_4S_2$
Purity	> 98%
MW	499.62
CAS Number	[130464-84-5]
SMILES	<chem>CC(C)CC(C(CSC1=CC=CS1)C(=O)N([O-]))C(=O)NC(CC2=CC=CC=C2)C(=O)NC.[Na+]</chem>
Expiration Date	12 months from date of receipt.