

Product Datasheet

Tipranavir (orb1182439)

Description	Tipranavir
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{31}H_{33}F_3N_2O_5S$
Purity	> 98%
MW	602.66
CAS Number	[174484-41-4]
SMILES	<chem>CCCC1(CC(=C(C(=O)O1)C(CC)C2=CC(=CC=C2)NS(=O)(=O)C3=NC=C(C=C3)C(F)(F)F)O)CCC4=CC=CC=C4</chem>
Expiration Date	12 months from date of receipt.