

## Product Datasheet

### Talnetant hydrochloride (orb1182442)

<b>Description</b>	Talnetant hydrochloride
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{25}H_{23}ClN_2O_2$
<b>Purity</b>	> 98%
<b>MW</b>	418.92
<b>CAS Number</b>	[204519-66-4]
<b>SMILES</b>	<chem>CCC(C1=CC=CC=C1)NC(=O)C2=C(C(=NC3=CC=CC=C32)C4=CC=CC=C4)O.Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.