

## Product Datasheet

### TAK-242 S enantiomer (orb1182443)

<b>Description</b>	TAK-242 S enantiomer
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{15}H_{17}ClFNO_4S$
<b>Purity</b>	> 98%
<b>MW</b>	361.82
<b>CAS Number</b>	[243984-10-3]
<b>SMILES</b>	<chem>CCOC(=O)C1=CCCCC1S(=O)(=O)NC2=C(C=C(C=C2)F)Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.