

Product Datasheet

SB-505124 HCl (orb1182449)

Description	SB-505124 HCl
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{20}H_{22}ClN_3O_2$
Purity	> 98%
MW	371.86
CAS Number	[356559-13-2]
SMILES	<chem>CC1=CC=CC(=N1)C2=C(N=C(N2)C(C)(C)C)C3=CC4=C(C=C3)OC(=O)Cl</chem>
Expiration Date	12 months from date of receipt.