

## Product Datasheet

### PF-543 Citrate (orb1182455)

<b>Description</b>	PF-543 Citrate
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{33}H_{39}NO_{11}S$
<b>Purity</b>	> 98%
<b>MW</b>	657.73
<b>CAS Number</b>	[1415562-83-2]
<b>SMILES</b>	<chem>CC1=CC(=CC(=C1)OCC2=CC=C(C=C2)CN3CCCC3CO)CS(=O)(=O)C4=CC=CC=C4.C(C(=O)O)C(CC(=O)O)(C(=O)O)O</chem>
<b>Expiration Date</b>	12 months from date of receipt.