

Product Datasheet

KN-92 (orb1182475)

Description	KN-92
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{24}H_{25}ClN_2O_3S$
Purity	> 98%
MW	456.98
CAS Number	[176708-42-2]
SMILES	<chem>CN(CC=CC1=CC=C(C=C1)Cl)CC2=CC=CC=C2NS(=O)(=O)C3=CC=C(C=C3)OC</chem>
Expiration Date	12 months from date of receipt.