


Product Datasheet

Hoechst 33258 analog 6 (orb1182491)

Description	Hoechst 33258 analog 6 is a analog of Hoechst stains(Hoechst 33258), which are part of a family of b...
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{33}H_{40}N_6O$
Purity	> 98%
MW	536.71
CAS Number	[129244-66-2]
SMILES	<chem>CC(C)(C)C1=CC(=C2NC3=C(N2)C=C(C=C3)C4=NC5=C(N4)C=C(C=C5)N6CCN(CC6)C=C(C1=O)C(C)(C)C</chem>
Expiration Date	12 months from date of receipt.



Hoechst
33258 analog
6 Chemical
Structur...