


Product Datasheet

Hoechst 33258 analog 2 (orb1182494)

Description	Hoechst 33258 analog 2
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₆ H ₂₆ N ₆
Purity	> 98%
MW	422.52
CAS Number	[23491-54-5]
SMILES	CC1=CC=CC(=C1)C2=NC3=C(N2)C=C(C=C3)C4=NC5=C(N4)C=C(C=C5)N6CCN(CC6)C
Expiration Date	12 months from date of receipt.



Hoechst
33258 analog
2 Chemical
Structur...