

Product Datasheet

CCT241533 hydrochloride (orb1182520)

Description	CCT241533 hydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{28}ClFN_4O_4$
Purity	> 98%
MW	478.94
CAS Number	[1431697-96-9]
SMILES	<chem>CC(C)(C1CNCC1NC2=NC(=C3C=C(C=CC3=O)F)NC4=CC(=C(C=C42)OC)OC)O.Cl</chem>
Expiration Date	12 months from date of receipt.