

Product Datasheet

Bifeprunox Mesylate (orb1182673)

Description	Bifeprunox Mesylate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{27}N_3O_5S$
Purity	> 98%
MW	481.56
CAS Number	[350992-13-1]
SMILES	<chem>CS(=O)(=O)O.C1CN(CCN1CC2=CC=CC(=C2)C3=CC=CC=C3)C4=CC=C C5=C4OC(=O)N5</chem>
Expiration Date	12 months from date of receipt.