

Product Datasheet

TPT-260 (Dihydrochloride) (orb1182706)

Description	TPT-260 (Dihydrochloride)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_8H_{14}Cl_2N_4S_3$
Purity	> 98%
MW	333.32
CAS Number	[2076-91-7]
SMILES	<chem>C1=C(SC(=C1)CSC(=N)N)CSC(=N)N.Cl.Cl</chem>
Expiration Date	12 months from date of receipt.