

## Product Datasheet

### Deoxyvasicine HCl (orb1182713)

<b>Description</b>	Deoxyvasicine HCl
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{11}H_{12}N_2 \cdot HCl$
<b>Purity</b>	> 98%
<b>MW</b>	208.69
<b>CAS Number</b>	[61939-05-7]
<b>SMILES</b>	<chem>C1CC2=NC3=CC=CC=C3CN2C1.Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.