

Product Datasheet

HA-1077 dihydrochloride (orb1182747)

Description	HA-1077 dihydrochloride
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{14}H_{17}N_3O_2S \cdot 2HCl$
Purity	> 98%
MW	364.29
CAS Number	[203911-27-7]
SMILES	<chem>C1CNCCN(C1)S(=O)(=O)C2=CC=CC3=C2C=CN=C3.Cl.Cl</chem>
Expiration Date	12 months from date of receipt.