

Product Datasheet

Flurbiprofen Axetil (orb1182784)

Description	Flurbiprofen Axetil
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{19}FO_4$
Purity	> 98%
MW	330.35
CAS Number	[91503-79-6]
SMILES	<chem>CC(C1=CC(=C(C=C1)C2=CC=CC=C2)F)C(=O)OC(C)OC(=O)C</chem>
Expiration Date	12 months from date of receipt.